



New Mask system for mice and rats

Simple to use

for mice and rats

mounted on surgery table

or placed on table

combi-vet®



New Mask System with inlet/outlet adaptor

Overview

Our new mask system simplifies the anesthesia in mice and rats. The separate gas inlet as well as the gas outlet shown on picture are divided and not together in a hose. That means you can turn the outlet tube for the gas scavenging in any position and it makes the surgery handling easy. There is no need of a diaphragms anymore and the gas will be evacuated from the inlet part.

You can use the mask system with the Mask stabilizer or it can be placed on our surgery table from our program.

1. STANDARD



Figure 1




Figure 2

2. MOUNTED ON OUR SURGERY TABLE



Figure 3



Figure 4

Order-Nr. PS-0310 Mask system for mice

Order-Nr. PS-0311 Mask system for rats